


Planmeca dental units

Perfect balance and reach for instruments




Our mid-bending balanced arms contribute to a short activation track for instruments – minimising stress on the shoulders



The mid-bending instrument arms remain outside the patient's field of vision




Each operator can fine-adjust the balance according to their preference



The asymmetric mounting of the instrument hose and balanced arm prevent the hose from getting stuck behind the arm




Sideway movements ensure perfect access to the oral cavity




Optimised balance lightens the instrument's weight – minimising strain



It is possible to release the hose for longer treatments – the balanced arm just needs to be pulled out a little bit to activate the hose



Detachable instrument arms enable easy cleaning and changing of the arms



Different units – the same balanced arms

The instrument console of all Planmeca dental units features ergonomic balanced instrument arms. The console can be equipped with up to five or six instruments, depending on the unit model. Except for the syringe, all other slots are interchangeable and can be equipped with any instrument from Planmeca's extensive selection.

All instrument hoses feature intelligent Plug-and-play quick connectors – allowing for a new instrument to be taken into use easily. When an instrument hose is attached to the console, the unit automatically selects the correct services and provides the correct default settings for that specific instrument. The over-the-patient delivery arm ensures light, precise, and well-balanced movements in all situations.



Planmeca
Compact™ i Classic



Planmeca Compact™ i Touch



Planmeca Sovereign® Classic

PLANMECA

Asentajankatu 6 | 00880 Helsinki | Finland | tel. +358 20 7795 500 | sales@planmeca.com | www.planmeca.com

Images may contain optional items not included in standard delivery. Available configurations and features may have country or area specific variations. Some products displayed above may not be available in all countries or areas. Rights for changes reserved.

Planmeca, All in one, Anamat Plus, Cobra, Comfy, DentreVac, Digital perfection, Economat Plus, Elegant, Flexy, Mini-dent, Perio Fresh, PlanEasyMill, Planmeca 4D, Planmeca AINO, Planmeca ARA, Planmeca CAD/CAM, Planmeca Chair, Planmeca Clarify, Planmeca Compact, Planmeca FIT, Planmeca Intra, Planmeca iRomexis, Planmeca Lumion, Planmeca Lumo, Planmeca Maximity, Planmeca Minea, Planmeca Minendo, Planmeca Minetto, Planmeca mRomexis, Planmeca Noma, Planmeca Online, Planmeca PlanCAD, Planmeca PlanCAM, Planmeca PlanClear, Planmeca PlanMill, Planmeca Planosil, Planmeca PlanPure, Planmeca PlanScan, Planmeca PlanView, Planmeca ProCeph, Planmeca ProFace, Planmeca ProID, Planmeca ProMax, Planmeca ProModel, Planmeca ProOne, Planmeca ProScanner, Planmeca ProSensor, Planmeca ProX, Planmeca Romexis, Planmeca Serenus, Planmeca SingLED, Planmeca Sovereign, Planmeca Ultra Low Dose, Planmeca Vision, Planmeca Verity, Planmeca Waterline Cleaning System, Planmeca Xtremity, Proline Dental Stool, ProTouch, Saddle Stool, SmartPan, SmartTouch, Trendy and Ultra Relax are registered or non-registered trademarks of Planmeca in various countries.

10040980/1016/en

