

Features

1 Three-way syringe (air/water/spray) with sterilisable tips

3 instrument hoses with light (6-hole connector)

Integrated spray water system for syringe and instruments (check valve prevents the back flow of water)

1 Spray water tank: 0.75 litres

Spray adjustment regulators for each instrument

1 High-performance suction hose 300 l/min. Continuously adjustable. Unique BPR Swiss technology
Adapter for saliva ejector tips and oral-dry-cup

1 Suction fluid container: 1 litre (automatic overflow protection)

Sturdy, air- and water-tight plastic construction, will not scratch, dent or corrode

Detachable front cover

Disc type foot control

Integrated compressor: Capacity 75 l/min at 4 bar;
Sound pressure level: 48 dB; Compressed air tank: 4 litres

Technical Data

Air consumption instruments: 35-55 l/min

Handpiece pressure: 2.5-3.5 bar

Each instrument is individually adjustable

Dimensions treatment position: 540mm x 530mm x 1070mm

Dimensions transport position: 540mm x 530mm x 760mm

Weight: 47 kg

Case unit dimensions: 370mm x 190mm x 500mm

Case unit weight: 8 kg

Warranty: 2 years

Options: Electric scaler, curing light, electric micromotor, air motor, turbine, contra-angles, arm with tray, arm with lamp, DENTA-TROLLEY, DENTA-CASE etc.

Combi-Port Premium Mobile comfort

The COMBI-PORT PREMIUM offers you maximum performance and optimal comfort for treatments outside your dental clinic.

All optional instruments can be operated at full efficiency in continuous use, such as micromotors, turbines, scalers or even air polishers. The high-performance COMBI-PORT PREMIUM is rounded off by its highly efficient aspirator of 300 l/min (unique BPR Swiss technology). The suction capacity is continuously adjustable.

Compactly fitted into a trolley, after work the COMBI-PORT PREMIUM can easily be collapsed and wheeled to the next workplace like a suitcase. The high-performance COMBI-PORT PREMIUM compressor not only enables working throughout the day without interruption, but is also distinguished by its remarkably quiet operation (48 dB!).

Self-contained high-performance suction 300 l/min

Full performance of all handpieces

Options Handpieces freely selectable

Silent oil-free high-performance compressor

-  Continuous Operation
-  Full Performance
-  Plug-and-Play
-  All-in-One
-  Mobile/Rollable



Option
30 freely selectable programs
incl. 10 Endo



Rolling easy transport